

Stealth Ultra Headrest Includes:

1. Occipital Pad
2. Sub Occipital Pad
3. Facial Lateral Pads
4. Link style adjustable channel
5. Headrest mounting hardware

TECHNICAL DATA

1. 5/32 tool for Link and Triangle adjustment.
2. 2x Links included; 1x In-Line, 1x Offset. Links measure 3" from screw to screw.
3. Headrest attaches to: HMO475, TWB485, TWB480
4. Set screw adjustment for increased stability.

SELLING POINTS

1. Fully adjustable fore and aft movement for infinite positioning on links.
2. Triangle has full rotation for offset.
3. Anti Rotation slot on mounting rod.
4. Swing away lateral pads or switches can be added for additional support or drive control.
5. Smooth finished edging on every piece.
6. Lever action handle for easy seat back mount adjustment.
7. Ball mounted for rotational adjustment.

BUILDING YOUR ULTRA

Step 1: Occipital Pads

Occipital pads come in five sizes with a variety of foam and cover options. According to the individual's needs, determine the necessary pad size, whether foam or gel would be beneficial, and cover material.

500—510—520—530—CPN350

Step 2: Sub-Occipital Pads

Sub-Occipital pads come in six sizes, each uniquely designed for specific needs, with a variety of foam and cover options.

600—610—620—630—634—637

Step 3: Swing-Away

Facial Lateral Swing-Aways for pad, switch, speaker sites, and chin control devices can be added to both the occipital and sub occipital pads. (See swing away options on back)

Step 4: Hardware

The final step is deciding what type of hardware to use for mounting the Ultra to your chair. (See hardware options on back)

FOAM/GEL


Open Cell Foam
Closed Cell Foam
Action Gel

COVERS
















Dartex
Reverse Dartex
Vinyl
Sheep Skin
Neoprene

HARDWARE

TWB480
TWB480-FDM
HMO475
TWB485



Pad Sizes

Sub-Occipital		Occipital		Facial Lateral Pads	
 8.5"x2.75" 600	 10.5"x2.75" 620	 6"x3.75" 500	 3"x3" 520	 Open cell foam 6" x 3.75" CPN350	 6.375" x 1.5" 930
 6.375" x 2.25" 610	 630 Small	 4.5"x2" 510	 2"x3" 530	 2"x3" 910	 6.375" x 2.25" 920
 634 Medium	 637 Large			 4.5"x2" 900	

CLINICAL GUIDELINES

The occipital pad cradles the head and is primarily a resting spot with nominal lateral support, although it also acts as a posterior boundary minimizing hyper-extension. The occipital pad is necessary for transportation as well as tilt & recline features.


Positioned at the base of the skull, and used in conjunction with the occipital pad, the sub-occipital pad provides additional lateral and posterior support for individuals that require mild to moderate assistance in maintaining midline positioning.







Facial lateral pads are preferable for individuals with considerable uncontrolled movement such as seizures, reflexes, or individuals with low tone or cognition. Facial pads can be used to facilitate proper head alignment and deter further deformity or shortening of muscles.

More information available at www.stealthproducts.com

Swing Away Options

SUS9-B
Swing Away Mount for Ultra Series Headrests (Required for mounting any facial lateral pads, joystick mounts, or switches).



 CN600 Micro Lite switch mount (Flex Tubing)	 SM600 Mini Proportional Joystick Mount	 SM630 ASL MEC Joy-
 F8765 Figure 8 Mount (allows for a second swing away rod to be added)	 FO650 Fiber Optic Mount	 CN600 Sip & Puff mount

Switches

MC210
Tash Mini Cup Switch







ML210
Tash Micro Lite Switch



Egg Switches

 ESYEL	 ESBLK
 ESRED	 ESBLU
 ESGRN	

Available Hardware

 TWB480	 TWB480_FDM
 TWB485	 HMO475