

New Product From Stealth



All Positioning System

The APS is the newest headrest at Stealth Products and offers a unique flexibility. The innovative memory metal allows an “adjust on demand” quality unlike most headrests. We offer several attaching options, making it complimentary to many of our head positioning systems. With five sizes available, all spectrums of shapes can be met. Combined with our TWB hardware, our APS is one of the most versatile systems on the market. When you need comfort and adjustability, depend on our All Positioning System.



AP380 Shown Includes:
Medium All Positioning Headrest (AP305)
TWB Hardware (TWB480)

Durable Pad is constructed of a layered Closed/Open Cell foam coated metal shell and covered with a soft, durable stretch material

Cover is removable and washable



All materials meet Cal.117 Flammability Standards

Quick Reference

Combination part numbers for the All Positioning System with Hardware

(Headrest mounts to hardware with a 17mm Ball)



	SIZE	PAD ONLY	HM0475	TWB480	TWB485	HM0400
Headrest	XSmall	AP150	AP175	AP180	AP185	AP100
	Small	AP250	AP275	AP280	AP285	AP200
	Medium	AP350	AP375	AP380	AP385	AP300
	Large	AP450	AP475	AP480	AP485	AP400
	XLarge	AP550	AP575	AP580	AP585	AP500