



+1 (512) 715-9995
 +1 (800) 965-9229 Toll Free
 +1 (512) 715-9954 Fax
 +1 (800) 806-1225 Toll Free Fax

104 John Kelly Drive,
 Burnet, TX 78611
 information@sw.stealthproducts.com
 www.stealthproducts.com
 facebook.com/StealthProducts
 twitter.com/stealthproducts

This order form uses Adobe Acrobat® interactive technology and we recommended you to use Acrobat Reader DC or later after downloading and saving the form to your local drive.

EPIc Seating Order Form

Fields or sections with asterisk (*) are required.

PO#:

Branch*: Date*:

Contact*:


E-mail*:

Phone*: Fax*:

Mark for: Order Type*: Quote Order
PO required if placing an order.

US Only. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement.

1.0 CLIENT SEAT AREA MEASUREMENTS*



- A. Hip Width: _____
- B. Total Seat Depth: _____
- C. Cushion to Iliac Crest: _____
- D. Buttock to Greater Trochanter: _____
- E. Cushion to Greater Trochanter: _____
- F. Iliac Crest to Axilla: _____
- G. Iliac Crest to top of Shoulder: _____

1.1 MOBILITY BASE INFORMATION*

Make/Model*:

Quickie®
Iris®

Leggero
Enzo™

ki Mobility
Focus CR

Invacare®
Solara® 3G

Frame Information*:

Frame Width*: Frame Depth*:

Note: Best when paired with foot-operated tilt-release option and non-reclining canes.

2.0 EPiC SEATING INFORMATION

Resistance Levels

Low Force (50N)

Low to Normal Tone-
Ideal for users with limited active movement who may also experience posterior pelvic tilt. Used to allow relaxation of postural muscles while providing an articulating posterior thoracic support or to allow for manual adjustment and accommodation of limited active movement or kyphosis without the use of tools.

High Force (100N)

Normal to High Tone-
Ideal for users with normal to highly active movement, such as dystonic users who may exhibit an extension pattern of spasticity. Provides an articulating posterior thoracic support which can absorb energy exerted by the increased tone.

Posterior Pelvic Support Sizes

Width Options

- Small:** Fits 10"-16" frame width
- Medium:** Fits 16"-21" frame width
- Large:** Fits 21"-26" frame width

Height:

The posterior pelvic support is 6 inches tall. The pad is adjustable along the cane, allowing the pad to rest between the iliac crest and the PSIS.

2.1 EPiC SEATING CONFIGURATION* K0108

Selections in this section **MUST** match information provided in Mobility Base information.

ASM: Assembly
REM ASM: Removable Assembly

Quickie Iris

Low Force

- EPiC Seating ASM, Quickie® Iris, Small, 50N..... \$1900
Part: ESAT090901050
- EPiC Seating ASM, Quickie® Iris, Medium, 50N..... \$1900
Part: ESAT090902050
- EPiC Seating ASM, Quickie® Iris, Large, 50N..... \$1900
Part: ESAT090903050

High Force

- EPiC Seating ASM, Quickie® Iris, Small, 100N..... \$1900
Part: ESAT090901100
- EPiC Seating ASM, Quickie® Iris, Medium, 100N..... \$1900
Part: ESAT090902100
- EPiC Seating ASM, Quickie® Iris, Large, 100N \$1900
Part: ESAT090903100

Leggero Enzo

Low Force

- EPiC Seating ASM, Leggero Enzo™, Small, 50N \$1900
Part: ESAT040401050
- EPiC Seating REM ASM, Leggero Enzo™, Small, 50N..... \$1900
Part: ESAT040404050
- EPiC Seating ASM, Leggero Enzo™, Medium, 50N \$1900
Part: ESAT040402050
- EPiC Seating REM ASM, Leggero Enzo™, Medium, 50N.... \$1900
Part: ESAT040405050
- EPiC Seating ASM, Leggero Enzo™, Large, 50N \$1900
Part: ESAT040403050
- EPiC Seating REM ASM, Leggero Enzo™, Large, 50N \$1900
Part: ESAT040406050

High Force

- EPiC Seating ASM, Leggero Enzo™, Small, 100N \$1900
Part: ESAT040401100
- EPiC Seating REM ASM, Leggero Enzo™, Small, 100N \$1900
Part: ESAT040404100
- EPiC Seating ASM, Leggero Enzo™, Medium, 100N \$1900
Part: ESAT040402100
- EPiC Seating REM ASM, Leggero Enzo™, Medium, 100N . \$1900
Part: ESAT040405100
- EPiC Seating ASM, Leggero Enzo™, Large, 100N \$1900
Part: ESAT040403100
- EPiC Seating REM ASM, Leggero Enzo™, Large, 100N \$1900
Part: ESAT040406100

Ki Focus CR

Low Force

- EPiC Seating ASM, Ki Focus CR, Small, 50N \$1900
Part: ESAT030501050
- EPiC Seating ASM, Ki Focus CR, Medium, 50N \$1900
Part: ESAT030502050
- EPiC Seating ASM, Ki Focus CR, Large, 50N \$1900
Part: ESAT030503050

High Force

- EPiC Seating ASM, Ki Focus CR, Small, 100N \$1900
Part: ESAT030501100
- EPiC Seating ASM, Ki Focus CR, Medium, 100N \$1900
Part: ESAT030502100
- EPiC Seating ASM, Ki Focus CR, Large, 100N \$1900
Part: ESAT030503100

Invacare Solara 3G

Low Force

- EPiC Seating ASM, Invacare® Solara 3G, Small, 50N..... \$1900
Part: ESAT021601050
- EPiC Seating ASM, Invacare® Solara 3G, Medium, 50N.... \$1900
Part: ESAT021602050
- EPiC Seating ASM, Invacare® Solara 3G, Large, 50N \$1900
Part: ESAT021603050

High Force

- EPiC Seating ASM, Invacare® Solara 3G, Small, 100N \$1900
Part: ESAT021601100
- EPiC Seating ASM, Invacare® Solara 3G, Medium, 100N . \$1900
Part: ESAT021602100
- EPiC Seating ASM, Invacare® Solara 3G, Large, 100N \$1900
Part: ESAT021603100

3.0 SEAT CUSHION SELECTION

Omit Cushion

For use with an aftermarket cushion

The **XXYY** portion of part numbers in this section directly relate to the dimensions of the cushion, where **XX** represents the width and **YY** represents the depth. Available sizes are shown on the corresponding cushion chart.

Skin Protection & Positioning

TRU-Comfort 2 SPP

- TRU-Comfort 2 SPP Cushion, 14"-21" wide, E2607 \$420
- TRU-Comfort 2 SPP Cushion, 22" wide and greater, E2608 \$525

Cushion Part: ST-TC2SPP**XXYY**
 Cushion Width: _____
 Cushion Depth: _____

Depth	Width									
	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"
14"										
15"										
16"										
17"										
18"										
19"										
20"										
21"										
22"										
24"										

Essence SPP

- Essence SPP Cushion, No Wedge, 10"-20" Wide, E2607 \$420
- Essence SPP Cushion, No Wedge, 22" Wide, E2608 \$525
- Essence SPP Cushion, Small Wedge, 10"-20" Wide, E2607 \$457
- Essence SPP Cushion, Small Wedge, 22" Wide, E2608 \$562
- Essence SPP Cushion, Medium Wedge, 10"-20" Wide, E2607..... \$488
- Essence SPP Cushion, Medium Wedge, 22" Wide, E2608..... \$593
- Essence SPP Cushion, Large Wedge, 10"-20" Wide, E2607..... \$520
- Essence SPP Cushion, Large Wedge, 22" Wide, E2608 \$625

Cushion Part: ST-ESSPP**XXYY**
 Cushion Width: _____
 Cushion Depth: _____
 Wedge: (Blank, S, M, or L) _____

Small wedge = 1/2" height
 Medium wedge = 1 1/2" height
 Large wedge = 2" height

Depth	Width										
	10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
10"											
12"											
14"											
15"											
16"											
17"											
18"											
19"											
20"											
21"											
22"											
24"											

Solution SPP

- Solution SPP Cushion, 10"-20" Wide, E2607..... \$394
- Solution SPP Cushion, 22"-24" Wide, E2608..... \$452

Cushion Part: ST-SOLSPP**XXYY**
 Cushion Width: _____
 Cushion Depth: _____

Depth	Width										
	10"	12"	14"	15"	16"	17"	18"	19"	20"	22"	24"
10"											
12"											
14"											
15"											
16"											
17"											
18"											
19"											
20"											
22"											
24"											

Spectrum Foam SPP

- Spectrum Foam SPP Cushion, 16"-20" Wide, E2607 \$483
- Spectrum Foam SPP Cushion, 22"-24" Wide, E2608 \$557
- Spectrum Foam SPP Cushion, 26"-32" Wide, E2608 \$620

Cushion Part: ST-SPCFSP**XXYY**
 Cushion Width: _____
 Cushion Depth: _____

Depth	Width											
	14"	16"	17"	18"	19"	20"	22"	24"	26"	28"	30"	32"
15"												
16"												
17"												
18"												
19"												
20"												
21"												
22"												
24"												
26"												
28"												

Spectrum Gel SPP

- Spectrum Gel SPP Cushion, 14"-20" Wide, E2607..... \$483
- Spectrum Gel SPP Cushion, 22"-24" Wide, E2608..... \$557

Cushion Part: ST-SPCGSP**XXYY**
 Cushion Width: _____
 Cushion Depth: _____

Depth	Width							
	14"	16"	17"	18"	19"	20"	22"	24"
15"								
16"								
17"								
18"								
19"								
20"								
21"								
22"								
24"								
26"								

Skin Protection

Glacial SP

- Glacial SP Cushion, 10"-20" Wide, E2603 \$347
- Glacial SP Cushion, 22"-24" Wide, E2604 \$462

Cushion Part: ST-GLSP**XXYY**
 Cushion Width: _____
 Cushion Depth: _____

Depth	Width											
	10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"
10"												
12"												
14"												
15"												
16"												
17"												
18"												
19"												
20"												
21"												
22"												
23"												
24"												

Seat Cushion Options continued on next page.

3.0 SEAT CUSHION SELECTION CONTINUED

ZEN SP

- ZEN SP Cushion, 10"-21" Wide, E2603 \$331
- ZEN SP Cushion, 22"-24" Wide, E2604 \$446

Cushion Part: ST-ZENSPXXYY
 Cushion Width: _____
 Cushion Depth: _____

		Width											
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"
Depth	10"												
	12"												
	14"												
	15"												
	16"												
	17"												
	18"												
	19"												
	20"												
	21"												
	22"												
	23"												
24"													

Positioning

Premiere P

- Premiere P Cushion, 14"-20" Wide, E2605 \$341
- Premiere P Cushion, 22" Wide, E2606 \$499

Cushion Part: ST-PREMPXXYY
 Cushion Width: _____
 Cushion Depth: _____

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Depth	13"								
	14"								
	15"								
	16"								
	17"								
	18"								
	19"								
	20"								
	22"								

General Use

Simplicity G

- Simplicity G Cushion, 10"-20" Wide, E2605 \$116
- Simplicity G Cushion, 22"-24" Wide, E2606 \$173

Cushion Part: ST-SIMGXXYY
 Cushion Width: _____
 Cushion Depth: _____

		Width											
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"
Depth	10"												
	12"												
	14"												
	15"												
	16"												
	17"												
	18"												
	19"												
	20"												
	21"												
	22"												
	23"												
24"													

4.0 BACK SUPPORT SELECTION*

Omit Backrest

For use with an aftermarket backrest

Choosing the right height back support: posterior support is divided into two segments, Pelvic and Thoracic. The 6" pelvic support is a standard component, please factor in a 2" gap between the two segments when selecting your thoracic back height.

Back Support Width*: _____

ADI Active (E2613)

- ADI AL Active Back Support, Low Height, 2Pt Pro Hardware (10") \$592.25
- ADI AL Active Back Support, Medium Height, 2Pt Pro Hardware (13") \$592.25
- ADI AL Active Back Support, Low Height, 4Pt Pro Hardware (10") \$592.25
- ADI AL Active Back Support, Medium Height, 4Pt Pro Hardware(13") \$592.25
- ADI AL Active Back Support, Low Height, Width Adjustable Elite Hardware (11") \$700.40
- ADI AL Active Back Support, Medium Height, Width Adjustable Elite Hardware (13") \$700.40

ADI Deep (E2615)

- ADI AL Deep Back Support, Low Height, 2Pt Pro Hardware (10") \$700.40
- ADI AL Deep Back Support, Medium Height, 2Pt Pro Hardware (13") \$700.40
- ADI AL Deep Back Support, Low Height, 4Pt Pro Hardware (10") \$700.40
- ADI AL Deep Back Support, Medium Height, 4Pt Pro Hardware (13") \$700.40
- ADI AL Deep Back Support, Low Height, Width Adjustable Elite Hardware (11") \$808.55
- ADI AL Deep Back Support, Medium Height, Width Adjustable Elite Hardware (13") \$808.55
- ADI AFT Deep Back Support, Low Height, 2Pt Pro Hardware (10") \$700.40
- ADI AFT Deep Back Support, Medium Height, 2Pt Pro Hardware (13") \$700.40
- ADI AFT Deep Back Support, Low Height, 4Pt Pro Hardware (10") \$700.40
- ADI AFT Deep Back Support, Medium Height, 4Pt Pro Hardware (13") \$700.40
- ADI AFT Deep Back Support, Low Height, Width Adjustable Elite Hardware (11") \$808.55
- ADI AFT Deep Back Support, Medium Height, Width Adjustable Elite Hardware (13") \$808.55

4.0 BACK SUPPORT SELECTION* CONTINUED

ADI Extra Deep (E2615)

ADI AL Extra Deep Back Support, Medium Height, 2Pt Pro Hardware (13")	\$700.40
ADI AL Extra Deep Back Support, Medium Height, 4Pt Pro Hardware (13")	\$700.40
ADI AL Extra Deep Back Support, Medium Height, Width Adjustable Elite Hardware (13")	\$808.55

5.0 CUSTOM BACK SELECTION (E2617)

Custom Back Selection

[See TFB Flex, TFB Contour or TFB Linear Order Forms.](#)

6.0 THORACIC LATERALS

Pad Options

Curved Lateral Pads (E0956)

UniLink 3x5 Lateral Pad, Neoprene, Metric, 1" Spacing	\$59.50
Part: ULP535050320	
UniLink 4x6 Lateral Pad, Neoprene, Metric, 1" Spacing	\$59.50
Part: ULP546050320	
UniLink 5x7 Lateral Pad, Neoprene, Metric, 1" Spacing	\$59.50
Part: ULP557050320	

Curved Lateral Pads, Gel (E0956)

UniLink 3x5 Gel Lateral Pad, Neoprene, Metric, 1" Spacing	\$98.00
Part: ULP535057120	
UniLink 4x6 Gel Lateral Pad, Neoprene, Metric, 1" Spacing	\$98.00
Part: ULP546057120	
UniLink 5x7 Gel Lateral Pad, Neoprene, Metric, 1" Spacing	\$98.00
Part: ULP557057120	

Flat Lateral Pads (E0956)

UniLink 3x5 Flat Lateral Pad, Neoprene, Metric, 1" Spacing	\$59.50
Part: ULP835050320	
UniLink 3x5 Flat Lateral Pad, CoolCore, Metric, 1" Spacing	\$59.50
Part: ULP835060320	
UniLink 3x5 Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing	\$59.50
Part: ULP835080320	
UniLink 4x6 Flat Lateral Pad, Neoprene, Metric, 1" Spacing	\$59.50
Part: ULP846050320	
UniLink 4x6 Flat Lateral Pad, CoolCore, Metric, 1" Spacing	\$59.50
Part: ULP846060320	
UniLink 4x6 Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing	\$59.50
Part: ULP846080320	
UniLink 5x7 Flat Lateral Pad, Neoprene, Metric, 1" Spacing	\$59.50
Part: ULP857050320	
UniLink 5x7 Flat Lateral Pad, CoolCore, Metric, 1" Spacing	\$59.50
Part: ULP857060320	
UniLink 5x7 Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing	\$59.50
Part: ULP857080320	

Flat Lateral Pads, Gel (E0956)

UniLink 3x5 Gel Flat Lateral Pad, Neoprene, Metric, 1" Spacing	\$98.00
Part: ULP835057120	
UniLink 3x5 Gel Flat Lateral Pad, CoolCore, Metric, 1" Spacing	\$98.00
Part: ULP835067120	
UniLink 3x5 Gel Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing	\$98.00
Part: ULP835087120	
UniLink 4x6 Gel Flat Lateral Pad, Neoprene, Metric, 1" Spacing	\$98.00
Part: ULP846057120	
UniLink 4x6 Gel Flat Lateral Pad, CoolCore, Metric, 1" Spacing	\$98.00
Part: ULP846067120	
UniLink 4x6 Gel Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing	\$98.00
Part: ULP846087120	

Thoracic Laterals continued on next page.

6.0 THORACIC LATERALS CONTINUED

Flat Lateral Pads, Gel (E0956)

UniLink 5x7 Gel Flat Lateral Pad, Neoprene, Metric, 1" Spacing	\$98.00
Part: ULP857057120	
UniLink 5x7 Gel Flat Lateral Pad, CoolCore, Metric, 1" Spacing	\$98.00
Part: ULP857067120	
UniLink 5x7 Gel Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing	\$98.00
Part: ULP857087120	

Thoracic Lateral Hardware

UniLink Hardware (E1028)

UniLink Swing-Away Lateral, Left Side	\$316.00	UniLink Swing-Away Lateral, Short, Left Side	\$316.00
Part: ULS25100-L		Part: ULS25000-L	
UniLink Swing-Away Lateral, Right Side	\$316.00	UniLink Swing-Away Lateral, Short, Right Side	\$316.00
Part: ULS25100-R		Part: ULS25000-R	

UniLink Mount Options (K0108)

No charge when ordering with the Lateral Hardware.

1 inch Cane Clamp	\$54.04
Part: UIMCNI100	
1 inch Crescent Clamp	\$54.04
Part: UIMCRI100	
90 Degree Lateral Bracket Mount (Planar, TFB I, Flex Mod.)	\$54.04
Part: UIMBKI090	
105 Degree Lateral Bracket Mount (TFB Rectangle)	\$54.04
Part: UIMBKI105	
135 Degree Lateral Bracket Mount (TFB Contour, TFB Flex Active, ADI, Q-Back)	\$54.04
Part: UIMBKI135	

Thin Series Hardware (E1028)

Lateral Thoracic Thin Swing-Away Hardware, Left Side	\$267.50
Part: TWSTL-L	
Lateral Thoracic Thin Swing-Away Hardware, Right Side	\$267.50
Part: TWSTL-R	
Lateral Thoracic Thin Swing-Away Hardware, Left Side (No Link)	\$267.50
Part: TWBTL-L	
Lateral Thoracic Thin Swing-Away Hardware, Right Side (No Link)	\$267.50
Part: TWBTL-R	

Thin Series Mount Options (K0108)

Lateral Thoracic Cane Mount, 1 inch	\$54.08
Part: LTHW500	

7.0 HEAD SUPPORTS

HEAD SUPPORT HARDWARE

Headrest Mounting Hardware (E1028)

Standard Straight Mounting Hardware, 17mm Ball, Metric	\$128.75
Part: HMO400-17M	
Multi-Axis Hardware, 17mm Ball, Metric	\$165.32
Part: HMO475-17M	
UniLink Headrest HW, 17mm Ball, 1/2" SQ 7" VB	\$296.50
Part: UHB2D33200 HCPCS: E1028	
UniLink Headrest HW, 17mm Ball, 15MM SQ 12" VB	\$296.50
Part: UHB3D33200 HCPCS: E1028	
UniLink Headrest HW, 17mm Ball, 1/2" SQ 12" VB	\$296.50
Part: UHB4D33200 HCPCS: E1028	

Headrest Hardware Mount Options (K0108)

Universal Headrest Adapter Plate	\$55.11
Part: SBM300U	

Head Support continued on next page.

7.0 HEAD SUPPORTS CONTINUED

HEADREST TONE DEFLECTOR

Headrest Tone Deflector, for 17MM Ball, Metric..... \$231.75
 Part: TD100M HCPCS: K0108

HEAD SUPPORT PADS

Comfort Plus (E0955)

Comfort Plus Nino 6" Headrest, Flat Plate, Metric..... \$170.25
 Part: CP06CCFPM

Comfort Plus Nino 8" Headrest, Flat Plate, Metric..... \$170.25
 Part: CP08CCFPM

Comfort Plus Nino 10" Headrest, Flat Plate, Metric..... \$170.25
 Part: CP10CCFPM

Comfort Plus Nino 14" Headrest, Flat Plate, Metric..... \$170.25
 Part: CP14CCFPM

Comfort Plus Nino 19" Headrest, Flat Plate, Metric..... \$259.96
 Part: CP19CCFPM

Adjustable Comfort Plus (E0955)

Adjustable Comfort Plus 14" Headrest, Metric..... \$338.34
 Part: ACP150M

Adjustable Comfort Plus 19" Headrest, Metric..... \$421.00
 Part: ACP550M

8.0 EPIC SEATING REPLACEMENT PARTS

Replacement Posterior Pelvic Block

Posterior Pelvic Block; 10 wide, 6in Tall..... \$145.00
 Part: ESR0008010

Posterior Pelvic Block; 11 wide, 6in Tall..... \$145.00
 Part: ESR0008011

Posterior Pelvic Block; 12 wide, 6in Tall..... \$145.00
 Part: ESR0008012

Posterior Pelvic Block; 13 wide, 6in Tall..... \$145.00
 Part: ESR0008013

Posterior Pelvic Block; 14 wide, 6in Tall..... \$145.00
 Part: ESR0008014

Posterior Pelvic Block; 15 wide, 6in Tall..... \$145.00
 Part: ESR0008015

Posterior Pelvic Block; 16 wide, 6in Tall..... \$145.00
 Part: ESR0008016

Posterior Pelvic Block; 17 wide, 6in Tall..... \$145.00
 Part: ESR0008017

Posterior Pelvic Block; 18 wide, 6in Tall..... \$145.00
 Part: ESR0008018

Posterior Pelvic Block; 19 wide, 6in Tall..... \$145.00
 Part: ESR0008019

Posterior Pelvic Block; 20 wide, 6in Tall..... \$145.00
 Part: ESR0008020

Posterior Pelvic Block; 21 wide, 6in Tall..... \$145.00
 Part: ESR0008021

Posterior Pelvic Block; 22 wide, 6in Tall..... \$145.00
 Part: ESR0008022

Posterior Pelvic Support

Lateral Pelvic Wedge Kit..... \$94.00
 Part: ESR0006000

Replacement Struts

Single Gas Strut, 50N \$154.00 ea.
 Part: ESR0003050

Single Gas Strut, 100N \$154.00 ea.
 Part: ESR0003100

Cable Assembly..... \$71.65 ea.
 Part: ESR0004000

Pivot Parts

Replacement Release Pin..... \$36.50 ea.
 Part: ESR0001000

Replacement Bearings \$82.50 ea.
 Part: ESR0002000

Attachment Tube

Invacare Solara 3G, Attachment Tube \$56.15 per pair.
 Part: ESR0200000

Ki Focus CR, Attachment Tube \$56.15 per pair.
 Part: ESR0300000

Leggero Enzo, Attachment Tube..... \$56.15 per pair.
 Part: ESR0400000

Quickie Iris, Attachment Tube \$56.15 per pair.
 Part: ESR0900000

Foam Grip Handles

Replacement Foam Grip Handles..... \$26.50 ea.
 Part: ESR0005000

NOTES
