
Stealth Tech Dashboard

SSP304 — 4 Point TF Mounting Hardware and Drop Stops w/ Solid Seat

The SSP304 Includes:

1. Solid Seat (Must specify all dimensions upon ordering.)
2. TF Mounts(x2)
3. Drop Stops(x2)



SSP304

Solid Seat with 4 Point Trouble Free Mounting Hardware and Drop Stops (E2609)

The SSP304 is a 1/2 inch wood solid seat with 4 Point Trouble Free Mounting Hardware and Drop Stops. The hardware is fully adjustable to mount to almost any chair size. The solid seat can be fully modified to accommodate the needs of almost any client. The cane clamps can be attached to 1", 7/8", or 3/4" canes (specify cane size)

Technical Data

1. 5/32 tool for 1/4– 20 button head screws (mounting hardware, adjustments).
2. 3/16 tool for 1/4—20 socket head screws (adjustments).


Selling Points

1. Easily removable knob & TF hardware makes it quick and simple to remove back.
2. Smooth finished edging on every piece.
3. Adjustable hardware to accommodate the needs of almost any client.


Standard Configuration Dimensions


Use these standard dimensions to determine which configuration you will need. See back for back modifications.

Cane clamp is designed to fit 1 inch, 7/8" inch or 3/4 inch canes. (Must specify upon ordering)



TF Mounts are also available in a smaller(space saver) option see seating hardware options section for part numbers & illustrations.






Mounting hardware has 3 inches of height adjustability to cane mount to almost any chair.


Stealth Tech Dashboard

SSP304 Modifications


Structural Seat Modifications




SMD315R
Right Side Offset




BMD223
"I" Seat



BMD223
"T" Seat



SMD320
Anti-Thrust Seat



SMD321
Contoured Seat


SMD317 – "I" Style Seat
SMD316 – "T" Style Seat
SMD318 – Built in Medial Thigh Support
SMD328 – Ischial Cutout for Seat
SMD357 – FIPS Package with Cover for Seat
SMD345 – Wedge Seat
SMD344 – Cross Brace Cut-out for Seat
SMD315L – Long Left Side Offset for Seat
SMD315R – Long Right Side Offset for Seat
SMD320 – Anti Thrust Addition for Seat
SMD321 – Contoured Seat
SMD336 – 3/8" Base for Seat (1/2" is Standard)
SMD362 – Flush Anti Thrust Addition for Seat
SMD379 – Cushion Separate from Wood Base


Seat Cover Modifications

SMD326 – Standard Cover for Seat
(Zipper Cover)
SMD360 – Velcro Cap Cover for Seat
SMD369 – Elastic Slip-On Cover for Seat
SMD309 – Monogram for Side of Seat

Custom Dimensions

All dimensions for solid seat and modifications must be specified upon ordering





Solid Seat

Foam Modifications

SMD300 – 1/2" Pudgee Foam
SMD301 – 1" Pudgee Foam
BMD204 – 1/2" Cubed Gel
SMD305 – 1/2" Extra Soft Sunmate
SMD306 – 1" Extra Soft Sunmate
SMD307– 1/2" Extra Soft T-Foam
SMD308 – 1" Extra Soft T-Foam
SMD310 – 1/2" Flat Gel
SMD311 – 1/2" Soft Sunmate
SMD312 – 1" Soft Sunmate
SMD313 – 1/2" Medium Sunmate
SMD314 – 1" Medium Sunmate
SMD330 – Grid Cut 2x2x2
SMD346 – 1/2" Soft T-Foam
SMD347 – 1" Soft T-Foam
SMD351 – 1/2" Medium T-Foam
SMD352 – 1" Medium T-Foam
SMD353 – 1/2" Firm Sunmate
SMD354 – 1" Firm Sunmate
SMD356 – 1/2" Cubed Gel
SMD359 – 1/4" Closed Cell Foam
SMD361 – 1/2" Closed Cell Foam
SMD363 – 1/2" HR1 Foam
SMD364 – 1" HR1 Foam
SMD365 – 1/2" HR2 Foam
SMD366 – 1" HR2 Foam
SMD367 – 1/2" HR3 Foam
SMD368 – 1" HR3 Foam